

SUMMIT MEDICAL

Equipment Company

UNIVERSAL CONTROL ARM (UCA) THE GOLD STANDARD FOR NON-REBREATHING SYSTEMS

PN: AS-01-0517



FEATURES:

1. Used with BAIN NRB System usually 60" to 72" long.
2. Reservoir Bag can be changed to match size of animal's tidal volume.
3. Independent "fine tune" Pressure Relief Valve which can be closed down to manually ventilate an animal via "bagging" by the technician, or mechanical IPPV (intermittent positive pressure ventilation).
4. Independent Manometer calibrated in cmH₂O (just like the rebreathing system). Safety: Using the manometer, the operator can see the approximate pressure in the animal's airway.
5. Concentration of anesthetic delivered can be changed almost instantaneously.
6. Can be used in special applications to safely anesthetize animals from 20gm. To 250 lb. Especially useful for MRI applications.



SUMMIT MEDICAL

Equipment Company

PRINCIPAL OF OPERATION:

1. Fresh gas from the O₂ flowmeter flows at a relatively high rate in relation to the animal's minute volume (tidal volume X breaths per minute).
2. Fresh gas flows directly to the distal end of the Bain Circuit that is attached to the animal's breathing apparatus either via endotracheal tube or face mask / nose cone.
3. During the apneic phase (the animal is not inhaling nor exhaling), fresh gas continues to flow at a relatively high rate. The fresh gas does not go into the animal's lungs, but flows down the exhaust side of the coaxial Bain system. This action flushes the last exhaled breath containing CO₂ down the exhaust side of the system. If the flowrate is sufficient, the area of the exhaust tube proximal to the animal is filled with fresh gas.
4. With the next inspiratory effort, the animal receives fresh gas from the fresh gas supply, and the volume of the exhaust side of the system which is now filled with fresh gas. If properly configured, the exhaust gas containing CO₂ does not reach the animal's airways again.
5. Using a long Bain NRB system (72"), it is possible to safely anesthetize an animal weighing 250 lbs. The limiting factor is the length of the exhaust side of the Bain and the flowrate for the oxygen flowmeter. In this case, the flowmeter would need to reach approximately 7-8 LPM.

CAUTION: DO NOT USE O₂ FLUSH THROUGH NRB SYSTEM

SUGGESTED O₂ FLOW RATES USING THE UNIVERSAL CONTROL ARM NON-REBREATHING SYSTEM (BAIN MANIFOLD WITH 72" ADULT BAIN NRB SYSTEM) BASED UPON ANIMAL WEIGHTS:

0 -	15 lbs	1.5 LPM O ₂
15 -	45 lbs	2.0 LPM O ₂
45 -	75 lbs	3.0 LPM O ₂
75 -	120 lbs	4.0 LPM O ₂
120 -	160 lbs	5.0 LPM O ₂
160 -	220 lbs	6.0 LPM O ₂
220 -	260 lbs	7.0 LPM O ₂

PLEASE NOTE: THESE ARE ONLY SUGGESTED STARTING POINTS FOR O₂ FLOW RATES. USE END TIDAL CO₂ MONITOR TO DETERMINE THE DEFINITIVE FLOW RATES.