

# SUMMIT MEDICAL

Equipment Company



Optional Wall Mount Platform



## **PAM** (Portable Anesthesia Machine)

Part Number: AS-01-0007 with the following:

- Portable, Lightweight, easy to operate
- Precision 0 - 4 LPM O2 flow meter with mechanical stop O2 flow control
- Silicone fresh gas tubing
- Vaporizer Inlet and Outlet adapter
- White powder coated aluminum
- Handle for carrying PAM
- Less than 1 cubic foot dimension
- Can be used with Universal Control Arm
- Available with new or reconditioned vaporizers (at extra cost)
- Available in Table Top or Wall Mount Platforms