

SUMMIT MEDICAL

Equipment Company

Posi-Vac Non-Rebreathing System

The Posi-Vac NRB System is used with in-house vacuum or any other vacuum source to completely remove waste anesthetic gases from the workplace and includes a Mapleson-D style body, nosecone with incisor bar and WAG collecting shroud. The Posi-Vac NRB System comes in three sizes and each nosecone/shroud sizes are interchangeable on the Mapleson-D style body. The Posi-Vac NRB Systems are compatible with any imaging modalities. Please call for more information.



X-Large Rat: AS-01-0510

Mouse/Rat Regular: AS-01-0512

Mouse: AS-01-0511